TalkLink Trust

	VOICEKEEDED	ACADELA	CDEAK LINIOLIE	MODELTALVER	CEDEVOICE	Annla vaicabanking
FOLUBA AFAIR ALTER TO		ACAPELA	•	MODELTALKER	CEREVOICE	Apple voicebanking
EQUIPMENT NEEDED	Laptop and headset OR directly onto Android and IOS Voicekeeper apps.	Laptop/computer and headset	Laptop/computer and headset	Laptop/computer and headset	Computer/laptop and headset	an iPhone, iPad or iMac device
Number of sentences to		50 minimum	150 minimum			
produce a voice Can add custom sentences to		350 maximum Yes – script of sentences	300 maximum No	298 Yes	100 No	150 No
improve quality and pronunciation	minimum 50	needs to be uploaded BEFORE you start recording.	NO	163	IVO	IVO
Ability to create voice from banked messages only	Yes	Yes	Yes – additional service called Voice Design. At least 5 minutes of speech recordings are needed to generate a successful voice	The audio quality must be very good in the recordings and there must be a sufficient amount recorded to be able to create a good quality voice	Yes, If good quality recordings, but this will cost extra	No
Time to complete voice banking	10 minutes – 1.5 hours	10 -60 minutes	1-3 hours	2 - 4 hours	30 minutes – 1 hour	30 minutes – 1 hour
Auditory prompts to support	Yes– sentence can be read	Yes - sentence can be	No - BUT they have a	Yes	Yes	No
literacy difficulty	out loud	read out loud	"simple word script" option providing easier, more accessible text to read. In addition, you can choose the script that you read, from varying levels of difficulty			
Ease of use	Easy and intuitive (website	1	Easy to follow the steps,	More timely and clunky	Not trialled by TalkLink	Quick and Easy
	preferred over phone application)	the steps	but takes longer to complete the process	to complete. Do not recommend if in a rush to complete.	clinicians due to cost	
Preview before purchase		Yes – can preview on website or Acapela app for 3 months trial	Yes – voice can be 'tweaked' upon request if user is not happy with it	compare three different versions of your synthetic	No	No purchase is required
Length of data storage	5 years	The voice created is stored for a period of one year in our system (legal commitment). On the practical side, this does not necessarily mean that we would delete the material, unless the patient or relatives instruct ACAPELA to do so.		voice Indefinitely (you consent to this when you sign up). people also can opt out of storage of their data; in which case their recordings and data are deleted within 30 days of their voice being delivered.	Indefinitely, but users need to pay before they sign up and start recording.	Stored directly to the persons iDevice
Cost (NZ)	apps and devices 2. approx. \$500 one off	 Ongoing subscription cost \$165 per year One off payment of \$1645 	1. \$350 for Voice build, 2. \$450 for Voice repair	\$165 NZ	\$1040	Free
Professionals given access to their own free voice, to show clients the comparisons	Yes	Yes	Yes	Yes	No	N/A
Compatible Communication	Large range of options.	Large range of options.	Limited options. Check	Check website for most	All Windows based	Not compatible with
apps and devices	Check website for most recent updates: https://www.thevoicekeeper.com/supported-aacdevices	Check website for most recent updates: https://mov.acapela-group.com/partners-products/products-apps/	website for most recent updates: https://preview.speakunique.co.uk/using-your-voice/	recent updates: <u>https://www.m</u>	programs (Grid 3; Communicator 5) Limited range of IOS apps	Windows or Android devices. IOS apps: Proloquo4text, Proloquo2Go, Pictello, Speak for Yourself,
Other Features	 Easy to do Website version preferred over phone app if possible Can choose to read 'short sentences' only if fatigue is an issue. Can make children's voices Very quick turn around User can record up to 2000 sentences if he wants to add in more custom script 'almost' hands free recording Voice can be installed as system voice on IOS 16 + Phone app works best with the latest IOS software. 	 Voice can be installed as system voice on IOS 16 + They also provide message banking support Can create a voice from your previously banked messages 	 Can do voice repair if speech quality has already declined App includes message banking service Regular online training sessions for professionals Can choose the genre and reading difficulty of the script you have to read 	service requiring user to wait for an email before progressing to the next step • ModelTalker is a project of the Nemours Speech Research Laboratory, part of a not-for-profit paediatric	voice on IOS 16 + Has not been tried or tested by any TalkLink Clinicians	Only works on iPhone 12 or later, iPad Air 5 th gen, iPad Pro 3 rd gen or later