

| | VOICEKEEPER | ACAPELA | SPEAK UNIQUE | MODELTALKER | CEREVOICE | Apple voicebanking |
|---|--|---|--|--------------------------------|---|-----------------------------------|
| EQUIPMENT NEEDED | Laptop and headset OR directly onto Android and IOS Voicekeeper apps. | Laptop/computer and headset | | Laptop/computer and headset | Computer/laptop and headset | an iPhone, iPad or iMac device |
| Number of sentences to produce a voice | 50 minimum 200 maximum | 50 minimum 350 maximum | 150 minimum 300 maximum | 298 | 100 | 150 |
| Can add custom sentences to improve quality and pronunciation | Yes – in addition to the minimum 50 | Yes – script of sentences needs to be uploaded BEFORE you start recording. | | | No | No |
| Ability to create voice from banked messages only | Yes | Yes | called Voice Design. At least 5 minutes of speech recordings are needed to generate a successful voice | recordings and there | Yes, If good quality recordings, but this will cost extra | No |
| Time to complete voice banking | 10 minutes – 1.5 hours | 10 -60 minutes | 1-3 hours | 2 - 4 hours | 30 minutes – 1 hour | 30 minutes – 1 hour |
| Auditory prompts to support literacy difficulty | Yes– sentence can be read out loud | Yes - sentence can be read out loud | No - BUT they have a "simple word script" option providing easier, more accessible text to read. In addition, you can choose the script that you read, from varying levels of difficulty | | Yes | No |

| | Easy and intuitive (website preferred over phone application) | | Easy to follow the steps, but takes longer to complete the process | More timely and clunky to complete. Do not recommend if in a rush to complete. | Not trialled by TalkLink clinicians due to cost | Quick and Easy |
|---|--|--|--|---|---|--|
| Preview before purchase | Yes - via website OR 30-day free trial on Voicekeeper IOS app | Yes – can preview on website or Acapela app for 3 months trial | Yes – voice can be 'tweaked' upon request if user is not happy with it | Yes –audition process where you get to compare three different versions of your synthetic voice | No | No purchase is required |
| Length of data storage | 5 years | The voice created is stored for a period of one year in our system (legal commitment). On the practical side, this does not necessarily mean that we would delete the material, unless the patient or relatives instruct ACAPELA to do so. | | Indefinitely (you consent to this when you sign up). people also can opt out of storage of their data; in which case their recordings and data are deleted within 30 days of their voice being delivered. | need to pay before they sign up and start | Stored directly to the persons iDevice |
| . , | 1.Free if using PRC- Saltillo apps and devices 2. approx. \$500 one off payment if done through the website 3. Recurring Monthly subscription via the app \$10 USD | cost \$165 per year 2. One off payment of | 1. \$350 for Voice build, 2. \$450 for Voice repair | \$165 NZ | \$1040 | Free |
| Professionals given access to their own free voice, to show clients the comparisons | Yes | Yes | Yes | Yes | No | N/A |
| Compatible Communication apps and devices | Large range of options. Check website for most recent updates: https://www.thevoicekee per.com/supported-aac- devices | recent updates: | Limited options. Check website for most recent updates: https://preview.speakuni que.co.uk/using-your-voice/ | recent updates: <u>https://www.m</u> | programs (Grid 3; Communicator 5) Limited range of IOS apps | Not compatible with Windows or Android devices. IOS apps: Proloquo4text, Proloquo2Go, Pictello, Speak for Yourself, |

| Other Features | • | Easy to do | • | Voice can be installed | • [| Can do voice repair if | • | More technical | • | Voice can be installed | Only works on |
|----------------|---|------------------------|---|------------------------|-----|------------------------|---|------------------------|---|------------------------|------------------------------------|
| | • | Website version | | as system voice on | | speech quality has | | service requiring user | ı | as system voice on | iPhone 12 or later, |
| | | preferred over phone | | IOS 16 + | | already declined | | to wait for an email | i | IOS 16 + | iPad Air 5 th gen, iPad |
| | | app if possible | • | They also provide | • | App includes | | before progressing to | • | Has not been tried or | Pro 3 rd gen or later |
| | • | Can choose to read | | message banking | | message banking | | the next step | 1 | tested by any | |
| | | 'short sentences' | | support | | service | • | ModelTalker is a | i | TalkLink Clinicians | |
| | | only if fatigue is an | • | Can create a voice | • | Regular online | | project of the | 1 | | |
| | | issue. | | from your previously | | training sessions for | | Nemours Speech | 1 | | |
| | • | Can make children's | | banked messages | | professionals | | Research Laboratory, | 1 | | |
| | | voices | | | • | Can choose the genre | | part of a not-for- | 1 | | |
| | • | Very quick turn | | | | and reading difficulty | | profit paediatric | 1 | | |
| | | around | | | | of the script you have | | healthcare system | 1 | | |
| | • | User can record up to | | | | to read | | | 1 | | |
| | | 2000 sentences if he | | | | | | | 1 | | |
| | | wants to add in more | | | | | | | 1 | | |
| | | custom script | | | | | | | 1 | | |
| | • | 'almost' hands free | | | | | | | 1 | | |
| | | recording | | | | | | | 1 | | |
| | • | Voice can be installed | | | | | | | 1 | | |
| | | as system voice on | | | | | | | 1 | | |
| | | IOS 16 + | | | | | | | i | | |
| | • | Phone app works | | | | | | | 1 | | |
| | | best with the latest | | | | | | | 1 | | |
| | | IOS software. | | | | | | | i | | |